mobee360

mobee metabolic

Montageanleitung
Installation manual

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info@mobee.eu • www.mobee.eu

1. mobee® metabolic device

The measuring device (mobee® metabolic hardware) is the METALYZER® 3B of the company CORTEX Biophysik GmbH from Leipzig.

The following information and notes are taken from the CORTEX product documentation and have been partially adapted for use in the context of resting metabolic rate measurement.

If you use the device in any other way than with mobee® metabolic the respective stricter requirements for maintenance, care, cleaning and disinfection apply!

If you have any questions, please contact your responsible distributor.

1.1 Contents

Overview contents:

Qty	Description
1x	METALYZER® 3B with mobee® front panel
1x	O2 sensor (installed in METALYZER® 3B)
1x	sample line, approx. 2m
1x	Flow sensor (turbine – disposable)
1x	CORTEX face mask set dark blue (S,M) incl. 2x adapters
1x	Nano Bluetooth USB adapter
1x	Head cap for face masks
1x	Medical desktop AC adapter
1x	Power cord Europe
1x	Gas-on-Demand-Ventil
1x	Disposable turbines PU 120
1x	Carrying case / box

1.2 About mobee® metabolic

The METALYZER® 3B is a spiroergometry device for stationary use. and works according to the Breath-by-Breath measurement method. In mobee® metabolic, only the functional range of the device required for resting metabolic rate measurement is used. During a resting metabolic rate measurement, the device transmits for each breath

- the oxygen concentration
- the carbon dioxide concentration
- as well as the volume flow

The values recorded by the device are transferred via USB interface to the PC and evaluated there using the mobee® metabolic software.

1.2. Device overview

No.	Connection	Funktion
1	gas	Connection for gas sample line
2	volume	Connection for flow sensor

No.	Display	LED	Function
3		blue	Display for Bluetooth® connection
4		green	Display during heart rate data reception
5		white	Display for operating status

No.	Switch	Function
6	on/off	on/off

No.	Connection	Function
7	Medical Power supply	Connection for medical desktop AC adapter
8	PC/USB	Connection USB cable

1.3 Start up

1.
Connect the device to a power outlet using the power supply unit/cable supplied. The connection may only be made to a power supply with a protective conductor.
Use only the supplied AC adapter and cable!

2.
Switch on the device by flipping the main switch. The white LED lights up. The device then performs a self-test, with the white LED flashing.
After switching on, the device requires at least 30 minutes to warm up until a stable operating temperature is reached and the measurement can be started!

3.
Connect the device to the PC. Use the supplied USB cable for this purpose. To do this, plug the wide connector of the USB cable into the „PC/USB“ socket on the back of the device. Then plug the flat connector of the USB cable into a free USB port of the computer.

Continuation to point 3.

Connect the volume flow sensor to the device by inserting the plug of the volume flow sensor into the socket marked „volume“.

Make sure that the color markings (red dot) on the plug and socket match.

To remove the plug, pull back the black outer cover of the plug a little and only then pull the plug carefully out of the socket.

Do not use force to lock or remove the plug, so as not to damage the pins or the plug itself!

4.

Connect the gas sample line to the device by inserting the pneumatic plug into the socket marked „gas“ on the device. Turn the plug clockwise until it clicks into place.

The plug will only fit in one direction. Proceed carefully and do not use force to avoid damaging the plug and socket!

**2. Assembling the face mask and connecting it to the mobee® metabolic device**

1.

Insert the adapter into the mask.

2.

Insert the disposable turbine into the adapter as far as it will go.

3.

Place the volume flow sensor on the disposable turbine. Turn the volume flow sensor to the right until it locks.

4.

Connect the sample line to the face mask by inserting the end of the sample line into the adapter.

5.

Connect the head cap with the face mask.

3. Inspection and calibration

In order to ensure valid measurements, the mobee® metabolic hardware must be calibrated regularly. This synchronizes the system with the ambient conditions.

A basic distinction is made between sensor testing using ambient air and sensor calibration using the calibration gas.

Caibration mode	Interval	Warm-up time
Sensor check with ambient air	After each device start or change of location	30-45 minutes
Sensor calibration with CORTEX calibration gas	1x monthly	30-45 minutes

3.1 Sensor check

Disconnect the gas aspiration section from the volume flow sensor, i.e. detach the gas sample line from the face mask and place the end of the gas sample line in a place where it is not exposed to breathing or ambient air. You can initiate the sensor check in the mobee® metabolic settings. Follow the instructions in the software.

Perform the sensor check each time you switch on the device or change the location of the device.

Please note that a valid sensor check is only possible with a sufficiently preheated device.

3.2 Sensor calibration

Sensor calibration is required **once a month** and is performed using the calibration gas. **Ensure sufficient warm-up time of the device before calibration.**

First prepare the calibration gas cyliner as described in the following chapter and then select the sensor calibration item in the mobee® metabolic settings. Follow the instructions in the software.

3.2.1 Preparing the calibration gas cyliner

1. Connect the gas-on-demand valve (A) with the adapter screwed on to the CORTEX calibration gas cylinder (B).

To do this, place the calibration gas cylinder on a flat surface and carefully screw the gas-on-demand valve straight onto the thread of the CORTEX calibration gas cylinder until you hear a short, quiet hiss.

Be careful not to tilt the gas-on-demand valve or damage the thread of the calibration gas cylinder.

2.

When used correctly, the display on the pressure valve of the gas-on-demand valve shows approx. 12 bar or just under 250 PSI.

3.

Insert the end of the gas sample line into the hose section on the gas-on-demand valve.

4. Support and contact

If you need help, please contact our support team.
(Prerequisite: valid service contract).

- Please submit support request via our [website](#)
- A ticket will be opened automatically in our support system.
- We will get back to you immediately by phone.

Support hours (CET)

Mon.–Fri. 9–12 a.m. and 1–4 p.m.
except lux. holidays & company closing days

Contact

SportMed SA – Headquarters
46, route de Wasserbillig
L-6490 Echternach

SportMed SA – Sales Office
91, Esplanade de la Moselle
L-6637 Wasserbillig

DE tel +49 6561 69771-72
LU tel +352 26720101

info@mobee.eu • www.mobee.eu

SportMed_{SA}

SportMed SA - 46, route de Wasserbillig
L-6490 Echternach - Luxembourg